


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Propyl acetate -  $\geq 99.5\%$

**Product Number:** 537438  
**CAS Number:** 109-60-4  
**MDL:** MFCD00009372  
**Formula:** C<sub>5</sub>H<sub>10</sub>O<sub>2</sub>  
**Formula Weight:** 102.13 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Refractive index at 25 ° C	1.380 - 1.384
Water (by Karl Fischer)	$\leq 0.1 \%$
Specific Gravity at 20 Degrees Celsius	0.885 - 0.890
Acidity as Acetic Acid	$\leq 0.01 \%$
Color Test, Pcs	$\leq 15$ APHA
Non-Volatile Residue g/100 mL	$\leq 0.01 \%$
Assay n-Propyl acetate	$\geq 99.50 \%$

Specification: PRD.0.ZQ5.10000080475

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.