

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

 β -Glucanase from *Trichoderma longibrachiatum*

Product Number:

G4423

CAS Number:

62213-14-3

TEST	Specification
------	---------------

% Protein (Biuret)

 ≥ 10

units/mg solid

 ≥ 1.0

Unit definition: One unit will liberate 1.0 umole of glucose from cellulose in one hr at pH 5.0 at 37°C (2 hr incubation time).

Specification: PRD.0.ZQ5.10000042762

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.