

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com Email USA: techsery@sial.com Outside USA: eurtechserv@sial.com

## **Product Specification**

Product Name:

Glycerol - for molecular biology, ≥99.0%

**Product Number:** G5516 CAS Number: 56-81-5

MDL: MFCD00004722

C3H8O3 Formula: Formula Weight: 92.09 g/mol

**TEST** Specification

Appearance (Turbidity) Clear Colorless Appearance (Color)

Appearance (Form) Liquid

Color Test < 10 APHA ≥ 99.0 % GC (area %)

Heavy Metal < 5 ppm as Lead

Solubility **Pass** 

Clear, colorless solution, 1 mL/mL in H2O. < 5 ppm

Elemental Analysis

Magnesium < 5 ppm Elemental Analysis

Protease by FITC-Casein None Detected DNAse, Exonuclease Detection None Detected RNAse Detection None Detected

NICKase, Endonuclease Detection None Detected Recommended Retest Period Confirmed

3 years

Specification: PRD.4.ZQ5.10000006262

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.