

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Urease from *Canavalia ensiformis* (Jack bean) - Type C-3, powder,  $\geq 600,000$  units/g solid

Product Number:

**U0251**

CAS Number:

9002-13-5

MDL:

MFCD00070858

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

units/G solid

 $\geq 600000$ 

Enzymatic Activity

Unit Definition: One micromolar unit will liberate 1.0 micromole of Ammonia from Urea per minute at pH 7.0 at 25 deg C.

Specification: PRD.0.ZQ5.10000056781

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.