

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

Silver carbonate on Celite® - extent of labeling: ~50 wt. % loading

Product Number:

363685



MDL:

MFCD00003403

Formula:

CAg2O3

Formula Weight:

275.75 g/mol

TEST**Specification**

Appearance (Color)

Faint Green to Green to Brown

Appearance (Form)

Conforms to Requirements

Powder and/or Chunks

Titration by KSCN

45.0 - 60.0 % wt

Specification: PRD.0.ZQ5.10000060944

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.