


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Maleimidoacetic acid N-hydroxysuccinimide ester - ≥95%

Product Number: M3884
CAS Number: 55750-61-3
MDL: MFCD00133523
Formula: C₁₀H₈N₂O₆
Formula Weight: 252.18 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear
10 mg/mL, CHCl ₃	
Proton NMR spectrum	Conforms to Structure
Purity (HPLC)	≥ 95 %

Specification: PRD.0.ZQ5.10000063083

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.