

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Methanol solution - contains 0.1 % (v/v) trifluoroacetic acid, 5 % (v/v) water, for HPLC

Product Number:

667390

CH₃OH

MDL:

MFCD00004595

Formula:

CH₄O

Formula Weight:

32.04 g/mol

TEST**Specification**

Appearance (Color)	Colorless
Appearance (Form)	Liquid
UV Absorbance 254nm	≤ 0.05
UV Absorbance 230nm	≤ 0.50
UV Absorbance 210nm	≤ 1.30
HPLC Absorbance 254 nm	≤ 10 mAU
Mass Balance	0.093 - 0.107 %
V% Trifluoroacetic Acid	
Mass Balance	4.5 - 5.5 %
V% Water	
Expiration Date Period	-----
18 months	

Specification: PRD.0.ZQ5.10000027219

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.