

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

17 α -Hydroxypregnenolone-20,21-¹³C₂,16,16-d₂ - 98 atom % D, 98 atom % ¹³C, 98% (CP)

Product Number:

803081

Formula:


13C₂C₁₉D₂H₃₀O₃

Formula Weight:

336.5 g/mol

Storage Temperature:

-20 °C



TEST

Specification

Appearance (Color)	White to Off White
Appearance (Form)	Powder
Purity (HPLC)	≥ 98 %
Proton NMR Spectrum	Conforms
Carbon NMR Spectrum	Conforms
Mass Spectrometry-D	≥ 98 %
Mass Spectrometry-13C	≥ 98 %

Specification: PRD.0.ZQ5.10000064923

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.