


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,3-Dimethyl-2-imidazolidinone - ≥99.0% (GC)

Product Number: 40727
CAS Number: 80-73-9
MDL: MFCD00003188
Formula: C₅H₁₀N₂O
Formula Weight: 114.15 g/mol



TEST

Specification

Appearance (Color)	Colorless to Very Faint Yellow
Appearance (Form)	Liquid
Infrared Spectrum	Conforms to Structure
Refractive index at 20 ° C	1.471 - 1.473
Water (by Karl Fischer)	≤ 0.2 %
Purity (GC)	≥ 99.0 %

Specification: PRD.0.ZQ5.10000050324

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.