

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Barium iodide dihydrate - ≥95%

Product Number: 223808
CAS Number: 7787-33-9
MDL: MFCD00149155
Formula: BaI₂ · 2H₂O
Formula Weight: 427.17 g/mol

BaI₂ · 2H₂O

TEST	Specification
------	---------------

Appearance (Color)	White to Tan
Appearance (Form)	Powder or Chunks
Titration by AgNO ₃ % I	≥ 56.4 %
ICP Major Analysis Confirms Barium Component	Confirmed

Specification: PRD.0.ZQ5.10000044104

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.