


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Iodo-5-methylthiophene - 97%

Product Number: 460532
CAS Number: 16494-36-3
MDL: MFCD00022492
Formula: C₅H₅IS
Formula Weight: 224.06 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)

Conforms to Requirements

Yellow to Yellow-Green to Dark Orange

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Assay

≥ 96.5 %

Specification: PRD.0.ZQ5.10000103071

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

