

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Cyanomethyl methyl(4-pyridyl)carbamodithioate - 98%

Product Number: **738689**
CAS Number: 1158958-89-4

Formula: C₉H₉N₃S₂
Formula Weight: 223.32 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
------	---------------

Appearance (Color)	White to Dark Yellow to Dark Orange
Appearance (Form)	Powder or Solid or Crystals or Crystalline Chunk(s) or Chunk(s)
Infrared spectrum	Conforms to Structure
Carbon	47.2 - 49.6 %
Nitrogen	18.3 - 19.3 %
Sulfur	28.0 - 29.4 %
Purity (HPLC)	≥ 97.5 %

Specification: PRD.2.ZQ5.10000032633

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.