

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-GABAA Receptor ($\alpha 3$ subunit) antibody produced in rabbit – affinity isolated antibody, lyophilized powder

Product Number:

G4291

MDL:

MFCD03455039

TEST	Specification
------	---------------

Formulation	Conforms
Lyophilized from phosphate buffered saline, pH 7.4, containing 1% BSA, 5% sucrose, and 0.025% sodium azide	
Protein Content	≥ 0.3
mg/mL after reconstitution	
Working Dilution by	Conforms
Direct Immunoblotting	
1:200 using anti-rabbit IgG peroxidase conjugate and ECL detection	
Immunohistochemistry	Conforms
1:200 using rat brain frozen sections	

Specification Date : 04/20/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.