


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Ethyl-1-propenyl ether, mixture of cis and trans - 98%

Product Number: 364045
CAS Number: 928-55-2
MDL: MFCD00015184
Formula: C₅H₁₀O
Formula Weight: 86.13 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
GC Purity (sum of isomers)	≥ 97.5 %A

Specification: PRD.0.ZQ5.10000026324

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.