

Product Specification

Product Name:

Monoclonal Anti-Synaptophysin antibody produced in mouse – clone SVP-38, ascites fluid

Product Number:

S5768

TEST	Specification
------	---------------

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Very Faint Yellow to Very Dark Yellow to Very Dark Brown
Protein by Biuret	15 - 55 mg/ml
Mouse IgG1 by RID	8.0 - 14.0 mg/ml
WB-Cell Line/Tissue Extract	Rat Brain
WB Titer	1:200-1:400
IHC - Tissue	Rat Cerebellum
IHC - Fixation	Frozen Sections or FFPE+ Antigen
IHC - Titer	Retrieval by Decloaking or Frozen Sections, Acetone Fixed
Recommended Retest Period	1:200-1:400
5 Years	-----

Specification: PRE.2.ZQ5.10000016148

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

