

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Isobutyric acid - ≥99%, FCC, FG

Product Number: W222208 CAS Number: 79-31-2

MDL: MFCD00002658

Formula: C4H8O2 Formula Weight: 88.11 g/mol

TEST

Appearance (Color)

Appearance (Form)
Refractive index at 20 ° C

Infrared spectrum

Purity (Titration by NAOH)

Specific Gravity

Reducing Substances Arsenic (As) Cadmium (Cd)

Mercury (Hg) Lead (Pb)

Expiration Date Period

5 Years

Specification: PRD.1.ZQ5.10000054423

Specification

Colorless

Liquid

1.392 - 1.395

Conforms to Structure

99.0 - 101.0 %

0.944 - 0.948

Pass

 \leq 3 ppm

< 1 ppm

 \leq 1 ppm

< 10 ppm

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.