


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Ethephon - ≥96% (titration)

Product Number: **C0143**
CAS Number: 16672-87-0
MDL: MFCD00002142
Formula: C₂H₆ClO₃P
Formula Weight: 144.49 g/mol
Storage Temperature: 2 - 8 °C



| TEST | Specification |
|------|---------------|
|------|---------------|

| | |
|----------------------------|-----------------------------|
| Appearance (Color) | White to Light Yellow |
| Appearance (Form) | Powder or Solid or Chunk(s) |
| Solubility (Color) | Colorless to Faint Yellow |
| Solubility (Turbidity) | Clear to Slightly Hazy |
| 25 mg/mL, H ₂ O | |
| Titration with NaOH | ≥ 96 % |

Specification: PRD.4.ZQ5.10000021055

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.