

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Rhodium on alumina - extent of labeling: 0.5 wt. % loading, pellets, 3.2 mm

Product Number:

206172

MDL:

MFCD00011201

Rh

TEST	Specification
------	---------------

Appearance (Color)

Conforms to Requirements

Grey

Appearance (Form)

Pellets

Rhodium

0.45 - 0.55 %

Loss on Drying

Specification: PRD.0.ZQ5.10000043580

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.