


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Propanol - ACS reagent, ≥99.5%

Product Number: 673773
CAS Number: 67-63-0
MDL: MFCD00011674
Formula: C₃H₈O
Formula Weight: 60.10 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 99.50 %
Acetone	≤ 0.002 %
Carbonyl compounds	≤ 0.002 %
as Propionaldehyde	
Water (by Karl Fischer)	≤ 0.2 %
Residue on Evaporation	≤ 0.001 %
Color Test	≤ 10 APHA
Solubility	Pass
Titrateable Acid or Base	≤ 0.0001 meq
Meets ACS Requirements	Current ACS Specification
Recommended Retest Period	-----
2 Years	

Specification: PRD.1.ZQ5.10000021097

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.