


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Bismuth(III) subsalicylate - 99.9% trace metals basis

Product Number: 480789
CAS Number: 14882-18-9
MDL: MFCD00085368
Formula: C₇H₅BiO₄
Formula Weight: 362.09 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Compleximetric EDTA % Bismuth	57.1 - 58.6 %
ICP Major Analysis Confirms Bismuth Component	Confirmed
Purity 99.9% Based On Trace Metals Analysis	Conforms
Trace Metal Analysis	≤ 1500.0 ppm

Specification: PRD.1.ZQ5.10000017169

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.