


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
5-Nitrofurfuryl alcohol - technical grade, 90%

Product Number: 433977
CAS Number: 2493-04-1
MDL: MFCD00272267
Formula: C₅H₅NO₄
Formula Weight: 143.10 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
Colorless to Brown	
Appearance (Form)	Liquid
Infrared Spectrum	Conforms to Structure
Purity (GC)	≥ 89.5 %

Specification: PRD.0.ZQ5.10000058962

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.