

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

2-AMINOETHYLISOTHIOURONIUM BROMIDE HYDRO

Product Number: A5879
CAS Number: 56-10-0
MDL: MFCD00037011
Formula: C₃H₉N₃S · 2HBr
Formula Weight: 281.01 g/mol

TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Purity (TLC)	≥ 98 %
Solubility (Turbidity) 50mg/ml in H ₂ O	Clear to Very Slightly Hazy
Solubility (Color)	Colorless to Yellow

Specification: PRD.0.ZQ5.10000039450

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.