


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Chlormadinone acetate - $\geq 98\%$

Product Number: **C5145**
CAS Number: 302-22-7
MDL: MFCD00056471
Formula: $C_{23}H_{29}ClO_4$
Formula Weight: 404.93 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Solubility (Color)
Solubility (Turbidity)
50 mg/mL, $CHCl_3$
Infrared spectrum
Purity (TLC)

White to Light Yellow
Powder
Colorless to Light Yellow
Clear

Conforms to Structure
 $\geq 98\%$

Specification: PRD.0.ZQ5.10000012218

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.