

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)


Email USA: [techserv@sial.com](mailto:techserv@sial.com)

Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Polyester-32-hydroxyl-1-acetylene bis-MPA dendron, generation 5 - 97%

Product Number: **686611**



Formula: C158H252O94  
Formula Weight: 3,653.6 AMU

TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Carbon	50.0 - 53.6 %
Molecular Weight	
Approximately 3654 g/mol	

Specification: PRD.0.ZQ5.10000035659

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.