


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Calcium phosphate monobasic - ≥95%

Product Number: **C8017**
CAS Number: 7758-23-8
MDL: MFCD00010898
Formula: $\text{H}_4\text{CaO}_8\text{P}_2$
Formula Weight: 234.05 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
50 mg/mL, 1 M HCl	
Loss on Drying	≤ 12 %
Titration by EDTA (Dry Basis)	≥ 95 %

Specification Date : 11/02/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.