

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Lead - powder, -100 mesh, 99.95% trace metals basis

Product Number: 391352
CAS Number: 7439-92-1
MDL: MFCD00134050
Formula: Pb
Formula Weight: 207.20 g/mol

Pb

TEST	Specification
Appearance (Color)	Conforms to Requirements
Grey to Black	
Appearance (Form)	Powder
Particle Size	Confirmed
-100 Mesh	
ICP Major Analysis	Confirmed
Confirms Lead Component	
Purity	Conforms
99.95% Based On Trace Metals Analysis	
Trace Metal Analysis	≤ 550.0 ppm

Specification: PRD.3.ZQ5.10000021109

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.