


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium phosphate monobasic - ReagentPlus®, ≥99.0%

Product Number: S0751
CAS Number: 7558-80-7
MDL: MFCD00003527
Formula: H₂NaO₄P
Formula Weight: 119.98 g/mol



| TEST | Specification |
|-----------------------------|--------------------|
| Appearance (Color) | White to Off White |
| Appearance (Form) | Powder or Crystals |
| Solubility (Color) | Colorless |
| Solubility (Turbidity) | Clear |
| 100 mg/ml, H ₂ O | |
| Loss on Drying | ≤ 1.0 % |
| Titration with NaOH | ≥ 99.0 % |
| Dry basis | |
| Recommended Retest Period | ----- |
| 4 years | |

Specification: PRD.2.ZQ5.10000045894

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.