


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**Product Name:
1-Pyrenemethanol 98%

Product Number: **389439**
CAS Number: 24463-15-8
MDL: MFCD00029252
Formula: C₁₇H₁₂O
Formula Weight: 232.28 g/mol

**TEST****Specification**

Appearance (Color) Off-White to Yellow to Yellow-Orange to Tan	Conforms to Requirements
Appearance (Form) Powder, Crystalline Powder, and/or Lumps	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Carbon	85.7- 90.1 %
Purity (HPLC)	≥97.5 %

Remarks:

Specification Date : 12/07/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.