

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Diisobutylaluminum hydride solution - 1.0 M in methylene chloride

Product Number:

214973

MDL:


MFCD00008928

Formula:

C₈H₁₉Al

Formula Weight:

142.22 g/mol

**TEST****Specification**

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Hydrogen Evolution

0.95 - 1.10 M

Specification: PRD.0.ZQ5.10000014785

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.