

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

Copper(II) perchlorate hexahydrate 98%

Product Number:**215392**

CAS Number:

10294-46-9

MDL:


MFCD00149666

Formula:

Cl₂CuO₈ · 6H₂O

Formula Weight:

370.54 g/mol



TEST	Specification
Appearance (Color)	Blue
Appearance (Form)	Conforms to Requirements
Compleximetric EDTA	16.7- 17.6 %
ICP Major Analysis	Confirmed

Remarks:

Specification Date : 21.06.2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.