


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Procainamide hydrochloride

Product Number: P9391
CAS Number: 614-39-1
MDL: MFCD00012998
Formula: C₁₃H₂₁N₃O · HCl
Formula Weight: 271.79 g/mol



TEST	Specification
Appearance (Color)	White to Yellow
Appearance (Form)	Powder
Titration by HClO ₄	≥ 98 %
Purity (TLC)	≥ 98 %
Solubility (Color)	Colorless to Yellow
Solubility (Turbidity)	Clear
50 mg/ml, H ₂ O	
Recommended Retest Period	-----
4 Years	

Specification: PRD.0.ZQ5.10000049022

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.