

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:


Potassium 2-fluoro-4-(trifluoromethyl)phenyltrifluoroborate - 90%

Product Number:

**731617**

CAS Number:

1150655-12-1



Formula:

C7H3BF7K

Formula Weight:

270.00 g/mol

**TEST****Specification**

Appearance (Color)

White to Beige

Appearance (Form)

Powder or Crystals or Chunk(s) or

Powder with Chunk(s) or Crystals with Chunk(s)

Fluorine NMR Spectra

Conforms to Structure

Boron HNMR Spectra

Conforms to Structure

Proton NMR spectra

Conforms to Structure

Carbon

27.40 - 34.88 %

Purity Titration NAOH by HCL

88.0 - 112.0 %

Specification: PRD.0.ZQ5.10000011543

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.