

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Lipase from porcine pancreas – Suitable for manufacturing of diagnostic kits and reagents, ≥20,000 units/mg protein, lyophilized powder

Product Number:

SRE0049

CAS Number:

9001-62-1

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

% Protein (Biuret)

40 - 70

units/mg protein

≥ 20000

(Using Olive oil) One unit will hydrolyze 1.0 micro-equivalent of fatty acid from a triglyceride in one hour at pH 7.7 at 37 deg C.

Recommended Retest Period

1 year

Specification: PRD.0.ZQ5.10000064739

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.