

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

NCTC 109 Medium - With sodium bicarbonate, without L-glutamine, liquid, sterile-filtered, suitable for cell culture

Product Number:

N1140

Storage Temperature:

2 - 8 °C

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Form)	Solution
pH	7.0 - 7.4
Osmolality	266 - 294 mOs/kg
Sterility by USP guidelines	Negative
Endotoxin Level	≤ 1.0 EU/ml
Glucose	0.9 - 1.9 g/l
Cell Culture Test	Pass
Cell Line	Attached Cell Line
L929 or 3T3	

Specification: PRD.0.ZQ5.10000046749

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.