

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

3-Chloro-4-pyridinecarboxaldehyde - 97%

Product Number:

636746

CAS Number:

72990-37-5

MDL:


MFCD06200712

Formula:

C₆H₄ClNO

Formula Weight:

141.56 g/mol



TEST

Specification

Appearance (Color)

White to Yellow

Appearance (Form)

Solid or Powder or Crystals or Chunks

Infrared spectrum

Conforms to Structure

Titration by Hydroxylamine

96.5 - 103.5 %

Purity (GC)

≥ 96.5 %

Specification: PRD.0.ZQ5.10000023371

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.