


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Clomipramine hydrochloride - $\geq 98\%$ (HPLC), powder

Product Number: C7291
CAS Number: 17321-77-6
MDL: MFCD00069234
Formula: C₁₉H₂₃ClN₂ · HCl
Formula Weight: 351.31 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear
25 mg/mL in H ₂ O	
Purity (HPLC)	$\geq 98\%$
Recommended Retest Period	-----
12 Years	

Specification: PRD.0.ZQ5.10000049903

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.