

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Insulin-Like Growth Factor-I antibody produced in goat – affinity isolated antibody

Product Number:

I8773

Storage Temperature:

-20 °C

TEST	Specification
Appearance (Form) Lyophilized from 250 µl of a 0.2 µm filtered solution in phosphate buffered saline pH 7.4 with 5% trehalose.	Lyophilized Powder
Bioactivity ND50 in presence of 6 ng/mL recombinant human IGF-1 using MCF-7 cells	2 - 12 µg/ml
Western Blot 0.1-0.2 µg/mL	Conforms
Endotoxin Assay < 0.1 EU/µg	Conforms

Specification: PRD.1.ZQ5.10000084355

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.