

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
JumpStart™ Taq ReadyMix™

Product Number: P2893

TEST	Specification
DNase, Exonuclease Detection	None Detected
RNAse Detection	None Detected
JumpStart Taq Antibody	Pass
Reduces non-specific product formation.	
Suitable for Quantitative PCR	Suitable

Specification: PRD.0.ZQ5.10000015446

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.