


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Agarose – For pulsed field electrophoresis running gel

Product Number: **A2929**
CAS Number: 9012-36-6
MDL: MFCD00081294



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Gel Temperature (10% gel)	38 - 43 °C
Gel Strength	≥ 2000 gm/cm ²
Electroendosmosis (-MR)	0.04 - 0.10
Sulfate (SO ₄)	≤ 0.20 %
Moisture Content	≤ 10 %
DNAse, Exonuclease Detection	None Detected
RNAse Detection	None Detected
NICKase, Endonuclease Detection	None Detected

Specification: PRD.3.ZQ5.10000009622

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.