

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Pectolyase from *Aspergillus japonicus* – lyophilized powder, ≥ 0.3 units/mg solid

Product Number:

P3026

MDL:

MFCD00131813

Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

Activity	≥ 0.3
Units/mg solid	
One unit will liberate 1.0 micromole galacturonic acid from polygalacturonic acid per min, pH 5.5, 25 deg C	

Specification: PRD.0.ZQ5.10000038839

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.