

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Ammonium molybdate tetrahydrate - BioReagent, suitable for cell culture, suitable for insect cell culture,

81.0-83.0% MoO3 basis

Product Number: M1019
CAS Number: 12054-85-2

CAS Number: 12054-85-2 MDL: MFCD00167059

Formula: H24Mo7N6O24 · 4H2O

Formula Weight: 1,235.86 g/mol

TEST Specification

Appearance (Color) White to Off White

Appearance (Form) Powder

Solubility (Color) Colorless to Very Faint Yellow Solubility (Turbidity) Clear

100 mg/mL, 1N NaOH

Infrared spectrum Conforms to Structure

Cell Culture Test Pass
Insect Cell Test Pass

Assay by Titration 81.0 - 83.0 %

as MoO3 by Pb(NO3)2

Recommended Retest Period 6 years

Specification: PRD.1.ZQ5.10000009973

(NH₄)₆Mo₇O₂₄ • 4H₂O

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.