

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Platinum(IV) oxide – surface area $\geq 60 \text{ m}^2/\text{g}$

Product Number:

206032



CAS Number:

1314-15-4

MDL:

MFCD00011184

Formula:

O₂Pt

Formula Weight:

227.08 g/mol

TEST

Specification

Appearance (Color)

Brown to Very Dark Brown

Appearance (Form)

Powder

ICP Major Analysis

Confirmed

Confirms Platinum Component

Gravimetric Analysis

 $\geq 80.0 \%$

% Platinum

Specification: PRD.0.ZQ5.10000028941

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.