


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium formate-d - 99 atom % D

Product Number: 373842
CAS Number: 3996-15-4
MDL: MFCD00075589
Formula: CDH₂NaO₂
Formula Weight: 71.03 g/mol



TEST	Specification
Water (by Karl Fischer)	≤ 0.50 %
IC	≥ 99.0 %
Proton NMR Spectrum	Conforms
Proton NMR	≥ 99.0 %
Deuterium NMR Spectrum	Conforms
¹³ C NMR Spectrum	Conforms
Appearance (Color)	White
Appearance (Form)	Powder

Specification: PRD.0.ZQ5.10000061081

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.