

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti-FLIP, N-Terminal antibody produced in rabbit - 1 mg/mL, affinity isolated antibody, buffered aqueous solution

**Product Number:** F9800

MDL: MFCD02686149

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Formulation	Conforms
-------------	----------

Solution in phosphate buffered saline containing 0.02% sodium azide.

Working Dilution	Conforms
------------------	----------

Immunoblotting: 0.5-1 ug/ml

Concentration	Conforms
---------------	----------

1 mg/ml

Specification: PRD.1.ZQ5.1000022171

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.