

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

1-Acetoxy-1,3-butadiene - mixture of cis and trans

**Product Number:****220868**

CAS Number:

1515-76-0

MDL:

MFCD00008714

Formula:


C<sub>6</sub>H<sub>8</sub>O<sub>2</sub>

Formula Weight:

112.13 g/mol

Storage Temperature:

2 - 8 °C

**TEST****Specification**

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Proton NMR spectra

Conforms to Structure

GC Purity (sum of isomers)

≥ 90.0 %A

Specification: PRD.1.ZQ5.10000034975

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.