


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Ethyl heptanoate – natural, ≥98%, FG

Product Number: W243728
CAS Number: 106-30-9
MDL: MFCD00009538
Formula: C₉H₁₈O₂
Formula Weight: 158.24 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.411 - 1.415
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 98.0 %
Specific Gravity	0.865 - 0.872
Acid Value	≤ 1.0 ml
Solubility (Turbidity)	Clear
Solubility (Color)	Colorless
1ml/3ml 70% Ethanol	
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm

Specification: PRD.2.ZQ5.10000017851

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.