

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Platinum(0)-1,3-divinyl-1,1,3,3-tetramethyldisiloxane complex solution - 0.05 M in poly(dimethylsiloxane), vinyl terminated

Product Number:

479535

MDL:


MFCD00151662

Formula:

C<sub>8</sub>H<sub>18</sub>O<sub>2</sub>PtSi<sub>2</sub>

Formula Weight:

381.48 g/mol



### TEST

### Specification

Appearance (Color)

Colorless to Yellow

Appearance (Form)

Viscous Liquid

Infrared spectrum

Conforms to Structure

ICP Major Analysis

Confirmed

Confirms Platinum Component

Assay

1.0 - 1.5 %

% Platinum by ICP

Specification: PRD.0.ZQ5.10000021347

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.