


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,3,5,6-Tetrafluoro-7,7,8,8-tetracyanoquinodimethane - 97%

Product Number: 376779
CAS Number: 29261-33-4
MDL: MFCD00042382
Formula: C₁₂F₄N₄
Formula Weight: 276.15 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
Yellow to Orange to Brown	
Appearance (Form)	Conforms to Requirements
Powder, Crystals, Crystalline Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
Carbon	50.4 - 54.0 %
Nitrogen	19.5 - 21.0 %
Purity (TLC)	≥ 97 %

Specification: PRD.0.ZQ5.10000015645

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.