

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Ketoconazole

Product Number:

UC280

CAS Number:

65277-42-1

Formula:


C₂₆H₂₈Cl₂N₄O₄

Formula Weight:

531.43 g/mol

Storage Temperature:

2 - 8 °C

**TEST****Specification**

Appearance (Color)

White to Light Yellow

Appearance (Form)

Solid

Purity (HPLC)

≥ 99.0 %

Proton NMR spectrum

Conforms to Structure

Specification: PRE.0.ZQ5.10000010795

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.