


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
tert-Butyl acetoacetate - reagent grade, 98%

Product Number: 537594
CAS Number: 1694-31-1
MDL: MFCD00008811
Formula: C₈H₁₄O₃
Formula Weight: 158.19 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Water (by Karl Fischer)	≤ 0.20 %
Color Test	≤ 100 APHA
Assay	≤ 2.00 %
T-Butyl Alcohol	
Miscellaneous Assay	≥ 98.00 %
Tert-Butyl Acetoacetate	

Specification: PRD.0.ZQ5.10000023933

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.